


# Siemens *syngo*® Media Viewer from Codonics

Siemens *syngo* Media Viewer allows the display of pre-aligned PET/CT or SPECT/CT images, as well as stand-alone CT, MR, SPECT or PET studies. Optimized for viewing fused studies, images are displayed in coronal, transaxial and sagittal planes with a correlated MIP, and fused images are displayed in a format that allows MIP blending between PET or SPECT images and the CT.



Codonics Virtua™ Medical Disc Publisher is the new standard for speed, efficiency and ease of use in automatic disc recorders. Virtua allows you to include the industry's leading viewers while its brilliant, full color disc labels increase your brand awareness.

For PET/CT

- View linked axial, coronal and sagittal views
- Navigate easily in three dimensions
- View MIP\* images correlated to coronal, transaxial and sagittal views
- Blend fused images
- Quantify the following: SUVs (standard uptake values), Hounsfield units, and voxel intensity

\* Maximum intensity projection

For more information contact Codonics at **800-444-1198** or **440-243-1198** and ask for our Virtua Sales Department.

Contact your Codonics representative for viewer configurations and availability. Specifications are subject to change without notice. Contact your Codonics representative for the latest information. Copyright © 2007 Codonics, Inc.



17991 Englewood Drive  
Middleburg Heights, OH 44130 USA  
(440) 243-1198  
(440) 243-1334 Fax  
Email [info@codonics.com](mailto:info@codonics.com)  
[www.codonics.com](http://www.codonics.com)